#### Abstract Submitted for the APR99 Meeting of The American Physical Society

Sorting Category: G.9 (Experimental)

Investigation of "White" Neutron Sources for Efficiency Calibration<sup>1</sup> T. N. MASSEY, A. SALAS BACCI, N. BOUKHAROUBA, C. E. BRIENT, S. M. CALHOUN, S. M. GRIMES, D. JACOBS, J. E. O'DONNELL, J. E. OLDENDICK, R. T. WHEELER, Ohio University, S. I. AL-QURASIHI, King Fahd University of Petroleum and Minerals, Saudi Arabia, W. E. PARKER, S.A. KREEK, K. RASCHKE, T. F. WANG, Lawrence Livermore National Laboratory — We have measured the neutrons spectrum produced by stopping targets of Al, B, BN, and Be with a deuteron beam using fission detectors and lithium glass detectors. The results of these measurements from 0.2 to 20 MeV will be presented. The limitations to accuracy of the determination of the neutron flux will be discussed. The advantages and disadvantages of targets to be used to produce a neutron spectrum to be used as a secondary standard will be discussed. The application to the measurement of the source reactions Be(p,n) and Be(d,n) will be demonstrated.

<sup>1</sup>Supported by DOE Grant DE-FG02-88ER40387



Prefer Oral Session Prefer Poster Session Thomas N. Massey massey@ohiou.edu Dept. Physic and Astr., Ohio University

Date submitted: December 4, 1998

Electronic form version 1.4

# Investigation of "White" Neutron Sources for Efficiency Calibration

T. N. Massey,

A. Salas Bacci, N. Boukharouba,
C. E. Brient, S. M. Calhoun,
S. M. Grimes, D. Jacobs,
J. E. O'Donnell, J. E. Oldendick,
R. T. Wheeler
Ohio University

S. I. Al-Quraishi King Fahd University of Petroleum and Minerals, Saudi Arabia

W. E. Parker, S.A. Kreek, K. Raschke, T. F. Wang Lawrence Livermore National Laboratory

#### Why Use a White Neutron Source for Calibration?

- Time Previous procedure took 1 week.
- Accuracy Previously used gas target reactions with at best 10% accuracy.
- New applied Physics, medical and neutron interogation, required better than 5% determination.
- This work is intended to generate an international secondary neutron standard.

#### Characteristic of an Ideal Source Reaction

- Lack of structure needed to avoid resolution dependent efficiencies.
- Need good neutron yield  $> 10^5$  neutrons/(MeV  $\mu$ C sr).
- A positive (d,n) Q-value is needed to cover the neutron energy range.
- The target should be stopping to avoid uncertainties due to target thickness.
- The target should be inexpensive, in high purity and stable in composition.
- Initial candidate targets were Be, B, BN, and Al.





BCI ADC 6

#### Data Reduction and Error Estimates

- The fission chamber with four thick  ${}^{235}$ U foils (  $0.9\% {}^{234}$ U,  $93.3\% {}^{235}$ U,  $0.3\% {}^{236}$ U, 5.5% ${}^{238}$ U ) was used for determination of the relative neutron flux and a thin  ${}^{235}$ U foil ( 1.0% ${}^{234}$ U,  $98.5\% {}^{235}$ U,  $0.43\% {}^{236}$ U,  $0.09\% {}^{238}$ U ) was used for absolute yield determination.
- ENDF/B-VI cross sections used to obtain average cross sections for correction factors and average fission cross sections. All cross sections used were convoluted with the resolution of the fission chamber.
- An additional effort was made to confirm the results of a previous mass spectroscopy result with an assay with gamma-rays. Contribution to the error in relative yield from the uncertainty in isotopic composition was less than 0.25%.
- The error (5%) in the absolute yield is dominated by the counting statistics of the thin standard <sup>235</sup>U foil.



## **Transmission Factors for Corrections**

















# Example Efficiency Determination Using Al(d, n)

- <sup>6</sup>Li Glass Detectors
- Stilbene Detector
- NE213 Detector





### Lithium Glass Detector Calibration





#### Conclusions

- We have measured neutron spectra from Be, B, BN and Al targets.
- We have used these results to determine efficiencies of neutron time-of-flight detectors.
- We have used at the beryllium, boron, and aluminum targets to obtain a calibration over the range of 40 keV to 20 MeV.
- Target Useful Energy Range
  - Be 0.050- $9.0~{\rm MeV}$ 
    - B 0.100- 19.0 MeV
  - BN None target deteriorates under irradiation
    - Al 0.200 14.0 MeV
      - Resultant efficiencies have been used to determine absolute neutron flux from the Be(p,n) and Be(d,n) reactions. The Be(p,n) reaction is currently being used by Harvard/MIT for Boron Neutron Capture Therapy (BNCT).